

# CD122/IL-2RB Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50451-1

## Capture Antibody Information

<b>Catalog Number:</b> 60174-3-PBS	<b>Clone ID:</b> 3F12D9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> interleukin 2 receptor, beta
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC025691	<b>Gene ID:</b> 3560
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag3929	

## Detection Antibody Information

<b>Catalog Number:</b> 60174-1-PBS	<b>Clone ID:</b> 3A12C7	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> interleukin 2 receptor, beta
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC025691	<b>Gene ID:</b> 3560
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag3929	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.098-25 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP50451-1 targets CD122/IL-2RB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD122/IL-2RB Monoclonal antibody, PBS Only (Capture/Detector) 60174-3-PBS (3F12D9). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD122/IL-2RB Monoclonal antibody, PBS Only (Detector) 60174-1-PBS (3A12C7). 100 µg. Concentration 1 mg/ml.

Alternative CD122/IL-2RB matched antibody pairs: MP01123-1, MP01123-2, MP01123-3, MP50451-2, MP50451-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50451-1, CD122/IL2RB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60174-3-PBS. Detection antibody: 60174-1-PBS. Standard: Ag3929. Range: 0.098-25 ng/mL.