

For Research Use Only

# FGFR1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50453-1

## Capture Antibody Information

Catalog Number:  
68982-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
2C3F7

Reactivity:  
human

GenBank:  
BC018128

Immunogen Catalog Number:  
Ag31801

Conjugate:  
Unconjugated

Full name:  
fibroblast growth factor receptor 1

Gene ID:  
2260

## Detection Antibody Information

Catalog Number:  
68982-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1B9C2

Reactivity:  
human

GenBank:  
BC018128

Immunogen Catalog Number:  
Ag31801

Conjugate:  
Unconjugated

Full name:  
fibroblast growth factor receptor 1

Gene ID:  
2260

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50453-1 targets FGFR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FGFR1 Monoclonal antibody, PBS Only (Capture/Detector) 68982-1-PBS (2C3F7). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: FGFR1 Monoclonal antibody, PBS Only (Detector) 68982-2-PBS (1B9C2). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50453-1, FGFR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68982-1-PBS. Detection antibody: 68982-2-PBS. Standard: Ag31801. Range: 0.781-100 ng/mL.