

For Research Use Only

Cyclin D1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50463-1

Capture Antibody Information

Catalog Number:
68986-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1G12F1

Reactivity:
human

GenBank:
BC000076

Immunogen Catalog Number:
Ag23228

Conjugate:
Unconjugated

Full name:
cyclin D1

Gene ID:
595

Detection Antibody Information

Catalog Number:
68986-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2G6A9

Reactivity:
human

GenBank:
BC000076

Immunogen Catalog Number:
Ag23228

Conjugate:
Unconjugated

Full name:
cyclin D1

Gene ID:
595

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50463-1 targets Cyclin D1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Cyclin D1 Monoclonal antibody, PBS Only (Capture) 68986-1-PBS (1G12F1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Cyclin D1 Monoclonal antibody, PBS Only (Detector) 68986-2-PBS (2G6A9). 100 μ g. Concentration 1 mg/ml.

Alternative Cyclin D1 matched antibody pairs: MP80003-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50463-1, Cyclin D1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68986-1-PBS. Detection antibody: 68986-2-PBS. Standard: Ag23228. Range: 0.098-100 ng/mL.