

SH2B3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50470-2

Capture Antibody Information

Catalog Number: 68989-3-PBS	Clone ID: 2C6A6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SH2B adaptor protein 3
Isotype: IgG1	GenBank: BC136451	Gene ID: 10019
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag17381	

Detection Antibody Information

Catalog Number: 68989-4-PBS	Clone ID: 1A1A7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SH2B adaptor protein 3
Isotype: IgG1	GenBank: BC136451	Gene ID: 10019
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag17381	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50470-2 targets SH2B3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SH2B3 Monoclonal antibody, PBS Only (Capture) 68989-3-PBS (2C6A6). 100 µg. Concentration 1 mg/ml.

Detection antibody: SH2B3 Monoclonal antibody, PBS Only (Detector) 68989-4-PBS (1A1A7). 100 µg. Concentration 1 mg/ml.

Alternative SH2B3 matched antibody pairs: MP50470-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50470-2, SH2B3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68989-3-PBS. Detection antibody: 68989-4-PBS. Standard: Ag17381. Range: 3.125-100 ng/mL.