

For Research Use Only

NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50472-2

Capture Antibody Information

Catalog Number:
68213-7-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A3A3

Reactivity:
human

GenBank:
BC064836

Immunogen Catalog Number:
Ag30868

Conjugate:
Unconjugated

Full name:
interleukin enhancer binding factor 3, 90kDa

Gene ID:
3609

Detection Antibody Information

Catalog Number:
68213-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3G6G6

Reactivity:
human

GenBank:
BC064836

Immunogen Catalog Number:
Ag30868

Conjugate:
Unconjugated

Full name:
interleukin enhancer binding factor 3, 90kDa

Gene ID:
3609

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50472-2 targets NF90,ILF3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NF90,ILF3 Monoclonal antibody, PBS Only (Capture) 68213-7-PBS (1A3A3). 100 µg. Concentration 1 mg/ml.

Detection antibody: NF90,ILF3 Monoclonal antibody, PBS Only (Detector) 68213-3-PBS (3G6G6). 100 µg. Concentration 1 mg/ml.

Alternative NF90,ILF3 matched antibody pairs: MP50472-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50472-2, NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68213-7-PBS. Detection antibody: 68213-3-PBS. Standard:Ag30868. Range: 0.781-50 ng/mL