

For Research Use Only

NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50472-3

Capture Antibody Information

Catalog Number:
68213-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2H10G3
Reactivity:
human, mouse, rat
GenBank:
BC064836
Immunogen Catalog Number:
Ag30868

Conjugate:
Unconjugated
Full name:
interleukin enhancer binding factor 3,
90kDa
Gene ID:
3609

Detection Antibody Information

Catalog Number:
68213-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1A9E3
Reactivity:
human
GenBank:
BC064836
Immunogen Catalog Number:
Ag30868

Conjugate:
Unconjugated
Full name:
interleukin enhancer binding factor 3,
90kDa
Gene ID:
3609

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50472-3 targets NF90,ILF3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NF90/ILF3 Monoclonal antibody, PBS Only (Capture) 68213-1-PBS (2H10G3). 100 µg. Concentration 1 mg/ml.

Detection antibody: NF90,ILF3 Monoclonal antibody, PBS Only (Detector) 68213-4-PBS (1A9E3). 100 µg. Concentration 1 mg/ml.

Alternative NF90,ILF3 matched antibody pairs: MP50472-1, MP50472-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50472-3, NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68213-1-PBS. Detection antibody: 68213-4-PBS. Standard:Ag30868. Range: 0.781-100 ng/mL