

For Research Use Only

CHCHD6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50477-3

Capture Antibody Information

Catalog Number:
66597-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2A10B5

Reactivity:
human

GenBank:
BC006123

Immunogen Catalog Number:
Ag14700

Conjugate:
Unconjugated

Full name:
coiled-coil-helix-coiled-coil-helix domain containing 6

Gene ID:
84303

Detection Antibody Information

Catalog Number:
66597-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B6D11

Reactivity:
human

GenBank:
BC006123

Immunogen Catalog Number:
Ag14700

Conjugate:
Unconjugated

Full name:
coiled-coil-helix-coiled-coil-helix domain containing 6

Gene ID:
84303

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50477-3 targets CHCHD6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHCHD6 Monoclonal antibody, PBS Only (Capture) 66597-2-PBS (2A10B5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CHCHD6 Monoclonal antibody, PBS Only (Detector) 66597-4-PBS (1B6D11). 100 μ g. Concentration 1 mg/ml.

Alternative CHCHD6 matched antibody pairs: MP50477-1, MP50477-2, MP50477-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50477-3, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-2-PBS. Detection antibody: 66597-4-PBS. Standard: Ag14700. Range: 0.781-100 ng/mL.