

For Research Use Only

CHCHD6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50477-4

Capture Antibody Information

Catalog Number:
66597-1-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
2A11E9
Reactivity:
human
GenBank:
BC006123
Immunogen Catalog Number:
Ag14700

Conjugate:
Unconjugated
Full name:
coiled-coil-helix-coiled-coil-helix domain containing 6
Gene ID:
84303

Detection Antibody Information

Catalog Number:
66597-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2G10C12
Reactivity:
human
GenBank:
BC006123
Immunogen Catalog Number:
Ag14700

Conjugate:
Unconjugated
Full name:
coiled-coil-helix-coiled-coil-helix domain containing 6
Gene ID:
84303

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50477-4 targets CHCHD6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHCHD6 Monoclonal antibody, PBS Only (Capture) 66597-1-PBS (2A11E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: CHCHD6 Monoclonal antibody, PBS Only (Detector) 66597-3-PBS (2G10C12). 100 µg. Concentration 1 mg/ml.

Alternative CHCHD6 matched antibody pairs: MP50477-1, MP50477-2, MP50477-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50477-4, CHCHD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66597-1-PBS. Detection antibody: 66597-3-PBS. Standard: Ag14700. Range: 0.391-100 ng/mL.