

For Research Use Only

TUFM Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50485-1

Capture Antibody Information

Catalog Number:
67802-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2C5B1

Reactivity:
human

GenBank:
BC010041

Immunogen Catalog Number:
Ag24976

Conjugate:
Unconjugated

Full name:
Tu translation elongation factor, mitochondrial

Gene ID:
7284

Detection Antibody Information

Catalog Number:
67802-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3F4D1

Reactivity:
human

GenBank:
BC010041

Immunogen Catalog Number:
Ag24976

Conjugate:
Unconjugated

Full name:
Tu translation elongation factor, mitochondrial

Gene ID:
7284

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50485-1 targets TUFM in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TUFM Monoclonal antibody, PBS Only (Capture) 67802-2-PBS (2C5B1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TUFM Monoclonal antibody, PBS Only (Detector) 67802-3-PBS (3F4D1). 100 μ g. Concentration 1 mg/mL.

Alternative TUFM matched antibody pairs: MP50485-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50485-1, TUFM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67802-2-PBS. Detection antibody: 67802-3-PBS. Standard: Ag24976. Range: 0.195-50 ng/mL