

For Research Use Only

FGFR2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50495-3

Capture Antibody Information

Catalog Number:
60403-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
3G11A2

Reactivity:
human

GenBank:

Conjugate:
Unconjugated

Full name:
fibroblast growth factor receptor 2

Gene ID:
2263

Detection Antibody Information

Catalog Number:
60403-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2C6G6

Reactivity:
human

GenBank:

Conjugate:
Unconjugated

Full name:
fibroblast growth factor receptor 2

Gene ID:
2263

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50495-3 targets FGFR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FGFR2 Monoclonal antibody, PBS Only (Capture) 60403-1-PBS (3G11A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: FGFR2 Monoclonal antibody, PBS Only (Detector) 60403-5-PBS (2C6G6). 100 µg. Concentration 1 mg/mL.

Alternative FGFR2 matched antibody pairs: MP01128-1, MP01128-2, MP01128-3, MP50495-1, MP50495-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50495-3, FGFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60403-1-PBS. Detection antibody: 60403-5-PBS. Standard:Eg0175. Range: 0.391-25 ng/mL