

For Research Use Only

Ninein Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50497-1

Capture Antibody Information

Catalog Number: 67132-2-PBS	Clone ID: 2G2E8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ninein (GSK3B interacting protein)
Isotype: IgG1	GenBank: BC034708	Gene ID: 51199
Purification Method: Protein G purification	Immunogen Catalog Number: Ag14652	

Detection Antibody Information

Catalog Number: 67132-3-PBS	Clone ID: 1E5A5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ninein (GSK3B interacting protein)
Isotype: IgG1	GenBank: BC034708	Gene ID: 51199
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag14652	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50497-1 targets Ninein in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Ninein Monoclonal antibody, PBS Only (Capture) 67132-2-PBS (2G2E8). 100 µg. Concentration 1 mg/mL.

Detection antibody: Ninein Monoclonal antibody, PBS Only (Detector) 67132-3-PBS (1E5A5). 100 µg. Concentration 1 mg/mL.

Alternative Ninein matched antibody pairs: MP50497-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50497-1, N9ein Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67132-2-PBS. Detection antibody: 67132-3-PBS. Standard: Ag14652. Range: 0.098-100 ng/mL.