

# CHAC1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50512-3

## Capture Antibody Information

<b>Catalog Number:</b> 60406-3-PBS	<b>Clone ID:</b> 2E2D9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> ChaC, cation transport regulator homolog 1 (E. coli)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC001847	<b>Gene ID:</b> 79094
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag7527	

## Detection Antibody Information

<b>Catalog Number:</b> 60406-4-PBS	<b>Clone ID:</b> 2E4H11	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> ChaC, cation transport regulator homolog 1 (E. coli)
<b>Isotype:</b> IgG2b	<b>GenBank:</b> BC001847	<b>Gene ID:</b> 79094
<b>Purification Method:</b> Protein A Magarose purification	<b>Immunogen Catalog Number:</b> Ag7527	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.195-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50512-3 targets CHAC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHAC1 Monoclonal antibody, PBS Only (Capture/Detector) 60406-3-PBS (2E2D9). 100 µg. Concentration 1 mg/ml.

Detection antibody: CHAC1 Monoclonal antibody, PBS Only (Detector) 60406-4-PBS (2E4H11). 100 µg. Concentration 1 mg/ml.

Alternative CHAC1 matched antibody pairs: MP50512-1, MP50512-2, MP50512-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50512-3, CHAC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60406-3-PBS. Detection antibody: 60406-4-PBS. Standard: Ag7527. Range: 0.195-100 ng/mL.