

IMMT Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50515-1

Capture Antibody Information

Catalog Number:
68226-1-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
2D2E8
Reactivity:
human, mouse, rat
GenBank:
BC002412
Immunogen Catalog Number:
Ag0102

Conjugate:
Unconjugated
Full name:
inner membrane protein,
mitochondrial (mitofilin)
Gene ID:
10989

Detection Antibody Information

Catalog Number:
68226-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2F10F2
Reactivity:
human
GenBank:
BC002412
Immunogen Catalog Number:
Ag0102

Conjugate:
Unconjugated
Full name:
inner membrane protein,
mitochondrial (mitofilin)
Gene ID:
10989

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50515-1 targets IMMT in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IMMT Monoclonal antibody, PBS Only (Capture/Detector) 68226-1-PBS (2D2E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: IMMT Monoclonal antibody, PBS Only (Detector) 68226-2-PBS (2F10F2). 100 µg. Concentration 1 mg/ml.

Alternative IMMT matched antibody pairs: MP50515-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50515-1, IMMT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68226-1-PBS. Detection antibody: 68226-2-PBS. Standard: Ag0102. Range: 0.098-100 ng/mL.