

For Research Use Only

SLC44A1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50519-2

Capture Antibody Information

Catalog Number:
60408-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1B10B6
Reactivity:
human
GenBank:
BC049203
Immunogen Catalog Number:
Ag6024

Conjugate:
Unconjugated
Full name:
solute carrier family 44, member 1
Gene ID:
23446

Detection Antibody Information

Catalog Number:
60408-4-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
1F2D5
Reactivity:
human
GenBank:
BC049203
Immunogen Catalog Number:
Ag6024

Conjugate:
Unconjugated
Full name:
solute carrier family 44, member 1
Gene ID:
23446

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50519-2 targets SLC44A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SLC44A1 Monoclonal antibody, PBS Only (Capture) 60408-3-PBS (1B10B6). 100 μ g. Concentration 1 mg/ml.
Detection antibody: SLC44A1 Monoclonal antibody, PBS Only (Detector) 60408-4-PBS (1F2D5). 100 μ g. Concentration 1 mg/ml.
Alternative SLC44A1 matched antibody pairs: MP50519-1
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50519-2, SLC44A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60408-3-PBS. Detection antibody: 60408-4-PBS. Standard: Ag6024. Range: 0.781-100 ng/mL.