

For Research Use Only

GPX4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50524-1

Capture Antibody Information

Catalog Number:
67763-2-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
1A1A5
Reactivity:
human
GenBank:
BC021567
Immunogen Catalog Number:
Ag30650

Conjugate:
Unconjugated
Full name:
glutathione peroxidase 4
(phospholipid hydroperoxidase)
Gene ID:
2879

Detection Antibody Information

Catalog Number:
67763-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2C3A1
Reactivity:
human
GenBank:
BC021567
Immunogen Catalog Number:
Ag30650

Conjugate:
Unconjugated
Full name:
glutathione peroxidase 4
(phospholipid hydroperoxidase)
Gene ID:
2879

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50524-1 targets GPX4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: GPX4 Monoclonal antibody, PBS Only (Capture) 67763-2-PBS (1A1A5). 100 µg. Concentration 1 mg/mL.
Detection antibody: GPX4 Monoclonal antibody, PBS Only (Detector) 67763-3-PBS (2C3A1). 100 µg. Concentration 1 mg/mL.
Alternative GPX4 matched antibody pairs: MP00049-1, MP00049-2, MP00049-3
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50524-1, GPX4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67763-2-PBS. Detection antibody: 67763-3-PBS. Standard: Ag30650. Range: 0.098-12.5 ng/mL