

For Research Use Only

SPRY2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50527-1

Capture Antibody Information

Catalog Number:
67123-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1E12C1
Reactivity:
human, mouse, rat
GenBank:
BC015745
Immunogen Catalog Number:
Ag28599

Conjugate:
Unconjugated
Full name:
sprouty homolog 2 (Drosophila)
Gene ID:
10253

Detection Antibody Information

Catalog Number:
67123-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1D9G1
Reactivity:
human
GenBank:
BC015745
Immunogen Catalog Number:
Ag28599

Conjugate:
Unconjugated
Full name:
sprouty homolog 2 (Drosophila)
Gene ID:
10253

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50527-1 targets SPRY2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SPRY2 Monoclonal antibody, PBS Only (Capture) 67123-1-PBS (1E12C1). 100 μ g. Concentration 1 mg/mL.
Detection antibody: SPRY2 Monoclonal antibody, PBS Only (Detector) 67123-2-PBS (1D9G1). 100 μ g. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50527-1, SPRY2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67123-1-PBS. Detection antibody: 67123-2-PBS. Standard: Ag28599. Range: 0.391-50 ng/mL