

For Research Use Only

SLC10A3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50537-1

Capture Antibody Information

Catalog Number:
60419-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H1A11

Reactivity:
human

GenBank:
BC004966

Immunogen Catalog Number:
Ag14038

Conjugate:
Unconjugated

Full name:
solute carrier family 10 (sodium/bile acid cotransporter family), member 3

Gene ID:
8273

Detection Antibody Information

Catalog Number:
60419-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
4A10B4

Reactivity:
human

GenBank:
BC004966

Immunogen Catalog Number:
Ag14038

Conjugate:
Unconjugated

Full name:
solute carrier family 10 (sodium/bile acid cotransporter family), member 3

Gene ID:
8273

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50537-1 targets SLC10A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC10A3 Monoclonal antibody, PBS Only (Capture) 60419-1-PBS (1H1A11). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC10A3 Monoclonal antibody, PBS Only (Detector) 60419-2-PBS (4A10B4). 100 µg. Concentration 1 mg/ml.

Alternative SLC10A3 matched antibody pairs: MP50537-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50537-1, SLC10A3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60419-1-PBS. Detection antibody: 60419-2-PBS. Standard: Ag14038. Range: 0.781-12.5 ng/mL