

For Research Use Only

LTF Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50548-4

Capture Antibody Information

Catalog Number:
68509-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
5G10E1

Reactivity:
human

GenBank:
BC015822

Immunogen Catalog Number:
Ag33645

Conjugate:
Unconjugated

Full name:
lactotransferrin

Gene ID:
4057

Detection Antibody Information

Catalog Number:
68509-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3D7B10

Reactivity:
human

GenBank:
BC015822

Immunogen Catalog Number:
Ag33645

Conjugate:
Unconjugated

Full name:
lactotransferrin

Gene ID:
4057

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50548-4 targets LTF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LTF Monoclonal antibody, PBS Only (Capture) 68509-2-PBS (5G10E1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: LTF Monoclonal antibody, PBS Only (Detector) 68509-3-PBS (3D7B10). 100 μ g. Concentration 1 mg/mL.

Alternative LTF matched antibody pairs: MP50548-5, MP50548-6

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50548-4, LTF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68509-2-PBS. Detection antibody: 68509-3-PBS. Standard: Ag33645. Range: 0.781-100 ng/mL.