

# LRPPRC Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50550-1

## Capture Antibody Information

<b>Catalog Number:</b> 67679-2-PBS	<b>Clone ID:</b> 3C8G4	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> leucine-rich PPR-motif containing
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC050311	<b>Gene ID:</b> 10128
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag30261	

## Detection Antibody Information

<b>Catalog Number:</b> 67679-3-PBS	<b>Clone ID:</b> 1F7B10	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> leucine-rich PPR-motif containing
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC050311	<b>Gene ID:</b> 10128
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag30261	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 3.125-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50550-1 targets LRPPRC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LRPPRC Monoclonal antibody, PBS Only (Capture/Detector) 67679-2-PBS (3C8G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: LRPPRC Monoclonal antibody, PBS Only (Detector) 67679-3-PBS (1F7B10). 100 µg. Concentration 1 mg/ml.

Alternative LRPPRC matched antibody pairs: MP50550-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50550-1, LRPPRC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67679-2-PBS. Detection antibody: 67679-3-PBS. Standard: Ag30261. Range: 3.125-100 ng/mL.