

For Research Use Only

CPVL Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50562-1

Capture Antibody Information

Catalog Number:
60429-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1D2G11

Reactivity:
human

GenBank:
BC016838

Immunogen Catalog Number:
Ag36547

Conjugate:
Unconjugated

Full name:
carboxypeptidase, vitellogenic-like

Gene ID:
54504

Detection Antibody Information

Catalog Number:
60429-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E6E8

Reactivity:
human

GenBank:
BC016838

Immunogen Catalog Number:
Ag36547

Conjugate:
Unconjugated

Full name:
carboxypeptidase, vitellogenic-like

Gene ID:
54504

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50562-1 targets CPVL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CPVL Monoclonal antibody, PBS Only (Capture) 60429-1-PBS (1D2G11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CPVL Monoclonal antibody, PBS Only (Detector) 60429-2-PBS (1E6E8). 100 μ g. Concentration 1 mg/mL.

Alternative CPVL matched antibody pairs: MP50562-2, MP50562-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50562-1, CPVL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60429-1-PBS. Detection antibody: 60429-2-PBS. Standard: Ag36547. Range: 3.125-100 ng/mL.