

For Research Use Only

PNPT1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50567-4

Capture Antibody Information

Catalog Number: 68309-5-PBS	Clone ID: 3F3D11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: polyribonucleotide nucleotidyltransferase 1
Isotype: IgG1	GenBank: BC053660	Gene ID: 87178
Purification Method: Protein G purification	Immunogen Catalog Number: Ag6290	

Detection Antibody Information

Catalog Number: 68309-3-PBS	Clone ID: 2E8H1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: polyribonucleotide nucleotidyltransferase 1
Isotype: IgG1	GenBank: BC053660	Gene ID: 87178
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6290	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-6.25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50567-4 targets PNPT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PNPT1 Monoclonal antibody, PBS Only (Capture) 68309-5-PBS (3F3D11). 100 µg. Concentration 1 mg/ml.

Detection antibody: PNPT1 Monoclonal antibody, PBS Only (Detector) 68309-3-PBS (2E8H1). 100 µg. Concentration 1 mg/ml.

Alternative PNPT1 matched antibody pairs: MP50567-1, MP50567-2, MP50567-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50567-4, PNPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68309-5-PBS. Detection antibody: 68309-3-PBS. Standard: Ag6290. Range: 0.098-6.25 ng/mL