

For Research Use Only

IL-12RB1/CD212 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50569-2

Capture Antibody Information

Catalog Number: 66829-2-PBS	Clone ID: 1F11F3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 12 receptor, beta 1
Isotype: IgG1	GenBank: BC029121	Gene ID: 3594
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6920	

Detection Antibody Information

Catalog Number: 66829-4-PBS	Clone ID: 1G2B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 12 receptor, beta 1
Isotype: IgG3	GenBank: BC029121	Gene ID: 3594
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag6920	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50569-2 targets IL-12RB1/CD212 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-12RB1/CD212 Monoclonal antibody, PBS Only (Capture) 66829-2-PBS (1F11F3). 100 µg. Concentration 1 mg/ml.

Detection antibody: IL-12RB1/CD212 Monoclonal antibody, PBS Only (Detector) 66829-4-PBS (1G2B5). 100 µg. Concentration 1 mg/ml.

Alternative IL-12RB1/CD212 matched antibody pairs: MP50569-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50569-2, IL-12RB1/CD212 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66829-2-PBS. Detection antibody: 66829-4-PBS. Standard: Ag6920. Range: 0.098-50 ng/mL.