

For Research Use Only

OMA1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50571-2

Capture Antibody Information

Catalog Number:
66572-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C4G4
Reactivity:
human
GenBank:
BC012915
Immunogen Catalog Number:
Ag11201

Conjugate:
Unconjugated
Full name:
OMA1 homolog, zinc metallopeptidase (*S. cerevisiae*)
Gene ID:
115209

Detection Antibody Information

Catalog Number:
67449-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E1D12
Reactivity:
human
GenBank:
BC012915
Immunogen Catalog Number:
Ag11201

Conjugate:
Unconjugated
Full name:
OMA1 homolog, zinc metallopeptidase (*S. cerevisiae*)
Gene ID:
115209

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50571-2 targets OMA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OMA1 Monoclonal antibody, PBS Only (Capture) 66572-3-PBS (2C4G4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: OMA1 Monoclonal antibody, PBS Only (Detector) 67449-2-PBS (2E1D12). 100 μ g. Concentration 1 mg/ml.

Alternative OMA1 matched antibody pairs: MP50571-1, MP50571-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50571-2, OMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66572-3-PBS. Detection antibody: 67449-2-PBS. Standard: Ag11201. Range: 0.781-100 ng/mL.