

For Research Use Only

Prx5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50572-2

Capture Antibody Information

Catalog Number:
67599-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1C1G10

Reactivity:
human

GenBank:
BC110983

Immunogen Catalog Number:
Ag12262

Conjugate:
Unconjugated

Full name:
peroxiredoxin 5

Gene ID:
25824

Detection Antibody Information

Catalog Number:
67599-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H1F2

Reactivity:
human

GenBank:
BC110983

Immunogen Catalog Number:
Ag12262

Conjugate:
Unconjugated

Full name:
peroxiredoxin 5

Gene ID:
25824

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50572-2 targets Prx5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Prx5 Monoclonal antibody, PBS Only (Capture/Detector) 67599-4-PBS (1C1G10). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Prx5 Monoclonal antibody, PBS Only (Detector) 67599-3-PBS (1H1F2). 100 μ g. Concentration 1 mg/ml.

Alternative Prx5 matched antibody pairs: MP50572-1, MP50572-3, MP50572-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50572-2, Prx5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67599-4-PBS. Detection antibody: 67599-3-PBS. Standard: Ag12262. Range: 0.098-100 ng/mL.