

ALCAM Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50574-2

Capture Antibody Information

Catalog Number:
67768-4-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
3H2G11

Reactivity:
human

GenBank:
BC137097

Immunogen Catalog Number:
Ag16892

Conjugate:
Unconjugated

Full name:
activated leukocyte cell adhesion molecule

Gene ID:
214

Detection Antibody Information

Catalog Number:
67768-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3C1D12

Reactivity:
human

GenBank:
BC137097

Immunogen Catalog Number:
Ag16892

Conjugate:
Unconjugated

Full name:
activated leukocyte cell adhesion molecule

Gene ID:
214

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50574-2 targets ALCAM in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ALCAM Monoclonal antibody, PBS Only (Capture) 67768-4-PBS (3H2G11). 100 µg. Concentration 1 mg/ml.

Detection antibody: ALCAM Monoclonal antibody, PBS Only (Capture/Detector) 67768-5-PBS (3C1D12). 100 µg. Concentration 1 mg/ml.

Alternative ALCAM matched antibody pairs: MP50574-1, MP50574-3, MP50574-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50574-2, ALCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67768-4-PBS. Detection antibody: 67768-5-PBS. Standard: Ag16892. Range: 6.25-200 ng/mL