

For Research Use Only

NGRN Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50577-1

Capture Antibody Information

Catalog Number:
68593-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2A10G9

Reactivity:
human

GenBank:
BC001682

Immunogen Catalog Number:
Ag6790

Conjugate:
Unconjugated

Full name:
neugrin, neurite outgrowth associated

Gene ID:
51335

Detection Antibody Information

Catalog Number:
68593-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F7C8

Reactivity:
human

GenBank:
BC001682

Immunogen Catalog Number:
Ag6790

Conjugate:
Unconjugated

Full name:
neugrin, neurite outgrowth associated

Gene ID:
51335

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50577-1 targets NGRN in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NGRN Monoclonal antibody, PBS Only (Capture) 68593-2-PBS (2A10G9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: NGRN Monoclonal antibody, PBS Only (Detector) 68593-3-PBS (1F7C8). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50577-1, NGRN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68593-2-PBS. Detection antibody: 68593-3-PBS. Standard: Ag6790. Range: 0.781-100 ng/mL.