

CYP11B1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50582-2

Capture Antibody Information

Catalog Number:
60435-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1C8A4

Reactivity:
human

GenBank:
BC096285

Immunogen Catalog Number:
Ag17158

Conjugate:
Unconjugated

Full name:
cytochrome P450, family 11,
subfamily B, polypeptide 1

Gene ID:
1584

Detection Antibody Information

Catalog Number:
60435-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2D12A2

Reactivity:
human

GenBank:
BC096285

Immunogen Catalog Number:
Ag17158

Conjugate:
Unconjugated

Full name:
cytochrome P450, family 11,
subfamily B, polypeptide 1

Gene ID:
1584

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50582-2 targets CYP11B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CYP11B1 Monoclonal antibody, PBS Only (Capture) 60435-3-PBS (1C8A4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CYP11B1 Monoclonal antibody, PBS Only (Detector) 60435-4-PBS (2D12A2). 100 µg. Concentration 1 mg/ml.

Alternative CYP11B1 matched antibody pairs: MP50582-1, MP50582-3, MP50582-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50582-2, CYP11B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60435-3-PBS. Detection antibody: 60435-4-PBS. Standard: Ag17158. Range: 0.098-50 ng/mL