

For Research Use Only

MRPS16 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50587-4

Capture Antibody Information

Catalog Number:
60439-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D12B12

Reactivity:
human

GenBank:
BC021106

Immunogen Catalog Number:
Ag9950

Conjugate:
Unconjugated

Full name:
mitochondrial ribosomal protein S16

Gene ID:
51021

Detection Antibody Information

Catalog Number:
60439-5-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A purification

Clone ID:
2A4A1

Reactivity:
human

GenBank:
BC021106

Immunogen Catalog Number:
Ag9950

Conjugate:
Unconjugated

Full name:
mitochondrial ribosomal protein S16

Gene ID:
51021

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50587-4 targets MRPS16 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MRPS16 Monoclonal antibody, PBS Only (Capture/Detector) 60439-3-PBS (1D12B12). 100 µg. Concentration 1 mg/ml.

Detection antibody: MRPS16 Monoclonal antibody, PBS Only (Detector) 60439-5-PBS (2A4A1). 100 µg. Concentration 1 mg/ml.

Alternative MRPS16 matched antibody pairs: MP50587-1, MP50587-2, MP50587-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50587-4, MRPS16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60439-3-PBS. Detection antibody: 60439-5-PBS. Standard: Ag9950. Range: 0.391-100 ng/mL.