

For Research Use Only

PLSCR3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50589-1

Capture Antibody Information

Catalog Number: 60440-1-PBS	Clone ID: 1B2F9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: phospholipid scramblase 3
Isotype: IgG1	GenBank: BC011735	Gene ID: 57048
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag13311	

Detection Antibody Information

Catalog Number: 60440-2-PBS	Clone ID: 1H3E12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: phospholipid scramblase 3
Isotype: IgG3	GenBank: BC011735	Gene ID: 57048
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag13311	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50589-1 targets PLSCR3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PLSCR3 Monoclonal antibody, PBS Only (Capture) 60440-1-PBS (1B2F9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: PLSCR3 Monoclonal antibody, PBS Only (Detector) 60440-2-PBS (1H3E12). 100 μ g. Concentration 1 mg/ml.

Alternative PLSCR3 matched antibody pairs: MP50589-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50589-1, PLSCR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60440-1-PBS. Detection antibody: 60440-2-PBS. Standard: Ag13311. Range: 0.195-100 ng/mL.