

For Research Use Only

SCOT Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50591-2

Capture Antibody Information

Catalog Number: 67836-2-PBS	Clone ID: 1F4F7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: 3-oxoacid CoA transferase 1
Isotype: IgG1	GenBank: BC009001	Gene ID: 5019
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag24792	

Detection Antibody Information

Catalog Number: 67836-4-PBS	Clone ID: 2E6C1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: 3-oxoacid CoA transferase 1
Isotype: IgG3	GenBank: BC009001	Gene ID: 5019
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24792	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50591-2 targets SCOT in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SCOT Monoclonal antibody, PBS Only (Capture) 67836-2-PBS (1F4F7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SCOT Monoclonal antibody, PBS Only (Detector) 67836-4-PBS (2E6C1). 100 μ g. Concentration 1 mg/mL.

Alternative SCOT matched antibody pairs: MP50591-1, MP50591-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50591-2, SCOT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67836-2-PBS. Detection antibody: 67836-4-PBS. Standard:Ag24792. Range: 0.391-100 ng/mL.