

MGME1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50593-2

Capture Antibody Information

Catalog Number:
67468-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F7A1

Reactivity:
human

GenBank:
BC016869

Immunogen Catalog Number:
Ag19644

Conjugate:
Unconjugated

Full name:
chromosome 20 open reading frame 72

Gene ID:
92667

Detection Antibody Information

Catalog Number:
67468-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1D7D1

Reactivity:
human

GenBank:
BC016869

Immunogen Catalog Number:
Ag19644

Conjugate:
Unconjugated

Full name:
chromosome 20 open reading frame 72

Gene ID:
92667

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50593-2 targets MGME1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MGME1 Monoclonal antibody, PBS Only (Capture) 67468-2-PBS (1F7A1). 100 µg. Concentration 1 mg/mL.

Detection antibody: MGME1 Monoclonal antibody, PBS Only (Detector) 67468-4-PBS (1D7D1). 100 µg. Concentration 1 mg/mL.

Alternative MGME1 matched antibody pairs: MP50593-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50593-2, MGME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67468-2-PBS. Detection antibody: 67468-4-PBS. Standard: Ag19644. Range: 6.25-200 ng/mL