

For Research Use Only

CNTROB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50605-3

Capture Antibody Information

Catalog Number:
67061-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3A8B9
Reactivity:
human
GenBank:
BC021134
Immunogen Catalog Number:
Ag25447

Conjugate:
Unconjugated
Full name:
centrobin, centrosomal BRCA2 interacting protein
Gene ID:
116840

Detection Antibody Information

Catalog Number:
67061-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3D11A1
Reactivity:
human
GenBank:
BC021134
Immunogen Catalog Number:
Ag25447

Conjugate:
Unconjugated
Full name:
centrobin, centrosomal BRCA2 interacting protein
Gene ID:
116840

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50605-3 targets CNTROB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: CNTROB Monoclonal antibody, PBS Only (Capture) 67061-5-PBS (3A8B9). 100 μ g. Concentration 1 mg/ml.
Detection antibody: CNTROB Monoclonal antibody, PBS Only (Capture/Detector) 67061-2-PBS (3D11A1). 100 μ g. Concentration 1 mg/ml.
Alternative CNTROB matched antibody pairs: MP50605-1, MP50605-2, MP50605-4
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50605-3, CNTROB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67061-5-PBS. Detection antibody: 67061-2-PBS. Standard:Ag25447. Range: 3.125-100 ng/mL.