

OCIAD1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50614-1

Capture Antibody Information

Catalog Number: 66698-1-PBS	Clone ID: 1C10C3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat	Full name: OCIA domain containing 1
Isotype: IgG2a	GenBank: BC003409	Gene ID: 54940
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9977	

Detection Antibody Information

Catalog Number: 66698-2-PBS	Clone ID: 2C3G12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: OCIA domain containing 1
Isotype: IgG1	GenBank: BC003409	Gene ID: 54940
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag9977	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50614-1 targets OCIAD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OCIAD1 Monoclonal antibody, PBS Only (Capture) 66698-1-PBS (1C10C3). 100 µg. Concentration 1 mg/ml.

Detection antibody: OCIAD1 Monoclonal antibody, PBS Only (Detector) 66698-2-PBS (2C3G12). 100 µg. Concentration 1 mg/ml.

Alternative OCIAD1 matched antibody pairs: MP50614-2, MP50614-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50614-1, OCIAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66698-1-PBS. Detection antibody: 66698-2-PBS. Standard: Ag9977. Range: 0.098-100 ng/mL.