

For Research Use Only

CHCHD5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50619-1

Capture Antibody Information

Catalog Number:
66584-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3F3B4

Reactivity:
human

GenBank:
BC004498

Immunogen Catalog Number:
Ag22534

Conjugate:
Unconjugated

Full name:
coiled-coil-helix-coiled-coil-helix domain containing 5

Gene ID:
84269

Detection Antibody Information

Catalog Number:
66584-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2H8B4

Reactivity:
human

GenBank:
BC004498

Immunogen Catalog Number:
Ag22534

Conjugate:
Unconjugated

Full name:
coiled-coil-helix-coiled-coil-helix domain containing 5

Gene ID:
84269

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50619-1 targets CHCHD5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHCHD5 Monoclonal antibody, PBS Only (Capture) 66584-2-PBS (3F3B4). 100 µg. Concentration 1 mg/mL.

Detection antibody: CHCHD5 Monoclonal antibody, PBS Only (Detector) 66584-3-PBS (2H8B4). 100 µg. Concentration 1 mg/mL.

Alternative CHCHD5 matched antibody pairs: MP50619-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50619-1, CHCHD5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66584-2-PBS. Detection antibody: 66584-3-PBS. Standard: Ag22534. Range: 0.098-50 ng/mL