

For Research Use Only

uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50621-2

Capture Antibody Information

Catalog Number:
60456-3-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
4C8B6

Reactivity:
human

GenBank:
NM_002659.4

Conjugate:
Unconjugated

Full name:
PLAUR

Gene ID:
5329

Detection Antibody Information

Catalog Number:
60456-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3C7D10

Reactivity:
human

GenBank:
NM_002659.4

Conjugate:
Unconjugated

Full name:
PLAUR

Gene ID:
5329

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50621-2 targets uPAR/CD87 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: uPAR/CD87 Monoclonal antibody, PBS Only (Capture) 60456-3-PBS (4C8B6). 100 µg. Concentration 1 mg/ml.

Detection antibody: uPAR/CD87 Monoclonal antibody, PBS Only (Detector) 60456-4-PBS (3C7D10). 100 µg. Concentration 1 mg/ml.

Alternative uPAR/CD87 matched antibody pairs: MP00129-1, MP00129-2, MP00129-3, MP00129-4, MP00129-5, MP50621-1, MP50621-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50621-2, uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60456-3-PBS. Detection antibody: 60456-4-PBS. Standard:Eg0901. Range: 0.098-25 ng/mL.