

# Porin Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50637-1

## Capture Antibody Information

**Catalog Number:**  
60463-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1F12D8  
**Reactivity:**  
human  
**GenBank:**  
BC008482  
**Immunogen Catalog Number:**  
Ag9241

**Conjugate:**  
Unconjugated  
**Full name:**  
voltage-dependent anion channel 1  
**Gene ID:**  
7416

## Detection Antibody Information

**Catalog Number:**  
60463-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2b  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1G7F5  
**Reactivity:**  
human  
**GenBank:**  
BC008482  
**Immunogen Catalog Number:**  
Ag9241

**Conjugate:**  
Unconjugated  
**Full name:**  
voltage-dependent anion channel 1  
**Gene ID:**  
7416

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50637-1 targets Porin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Porin Monoclonal antibody, PBS Only (Capture) 60463-1-PBS (1F12D8). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: Porin Monoclonal antibody, PBS Only (Detector) 60463-2-PBS (1G7F5). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative Porin matched antibody pairs: MP50637-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50637-1, Porin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60463-1-PBS. Detection antibody: 60463-2-PBS. Standard: Ag9241. Range: 0.098-100 ng/mL.