

Integrin beta-5 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50641-3

Capture Antibody Information

Catalog Number:
60471-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A10C7

Reactivity:
human

GenBank:
BC006541

Immunogen Catalog Number:
Ag29817

Conjugate:
Unconjugated

Full name:
integrin, beta 5

Gene ID:
3693

Detection Antibody Information

Catalog Number:
60471-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H1B2

Reactivity:
human

GenBank:
BC006541

Immunogen Catalog Number:
Ag29817

Conjugate:
Unconjugated

Full name:
integrin, beta 5

Gene ID:
3693

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50641-3 targets Integrin beta-5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Integrin beta-5 Monoclonal antibody, PBS Only (Capture) 60471-1-PBS (1A10C7). 100 µg. Concentration 1 mg/ml.

Detection antibody: Integrin beta-5 Monoclonal antibody, PBS Only (Detector) 60471-4-PBS (1H1B2). 100 µg. Concentration 1 mg/ml.

Alternative Integrin beta-5 matched antibody pairs: MP00512-1, MP00512-2, MP00512-3, MP00512-4, MP50641-1, MP50641-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50641-3, Integrin beta-5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60471-1-PBS. Detection antibody: 60471-4-PBS. Standard: Ag29817. Range: 0.781-100 ng/mL.