

For Research Use Only

HSD17B10 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50654-1

Capture Antibody Information

Catalog Number: 60069-2-PBS	Clone ID: 1C2C4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: hydroxysteroid (17-beta) dehydrogenase 10
Isotype: IgG1	GenBank: BC008708	Gene ID: 3028
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag1020	

Detection Antibody Information

Catalog Number: 60069-3-PBS	Clone ID: 1C9F2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: hydroxysteroid (17-beta) dehydrogenase 10
Isotype: IgG1	GenBank: BC008708	Gene ID: 3028
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag1020	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50654-1 targets HSD17B10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HSD17B10 Monoclonal antibody, PBS Only (Capture) 60069-2-PBS (1C2C4). 100 µg. Concentration 1 mg/ml.

Detection antibody: HSD17B10 Monoclonal antibody, PBS Only (Detector) 60069-3-PBS (1C9F2). 100 µg. Concentration 1 mg/ml.

Alternative HSD17B10 matched antibody pairs: MP50654-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50654-1, HSD17B10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60069-2-PBS. Detection antibody: 60069-3-PBS. Standard: Ag1020. Range: 0.098-100 ng/mL