

For Research Use Only

HSD17B10 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50654-2

Capture Antibody Information

Catalog Number:
60069-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1G5E5

Reactivity:
human, mouse, rat

GenBank:
BC008708

Immunogen Catalog Number:
Ag1020

Conjugate:
Unconjugated

Full name:
hydroxysteroid (17-beta)
dehydrogenase 10

Gene ID:
3028

Detection Antibody Information

Catalog Number:
60069-4-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
5F2F6

Reactivity:
human

GenBank:
BC008708

Immunogen Catalog Number:
Ag1020

Conjugate:
Unconjugated

Full name:
hydroxysteroid (17-beta)
dehydrogenase 10

Gene ID:
3028

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50654-2 targets HSD17B10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ERAB Monoclonal antibody, PBS Only (Capture) 60069-1-PBS (1G5E5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: HSD17B10 Monoclonal antibody, PBS Only (Detector) 60069-4-PBS (5F2F6). 100 μ g. Concentration 1 mg/ml.

Alternative HSD17B10 matched antibody pairs: MP50654-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50654-2, HSD17B10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60069-1-PBS. Detection antibody: 60069-4-PBS. Standard: Ag1020. Range: 0.098-50 ng/mL.