

For Research Use Only

PRAS40 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50656-1

Capture Antibody Information

Catalog Number:
60477-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F2A9

Reactivity:
human

GenBank:
BC051844

Conjugate:
Unconjugated

Full name:
AKT1 substrate 1 (proline-rich)

Gene ID:
84335

Detection Antibody Information

Catalog Number:
60478-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
1F11C10

Reactivity:
human

GenBank:
BC051844

Immunogen Catalog Number:
Ag15355

Conjugate:
Unconjugated

Full name:
AKT1 substrate 1 (proline-rich)

Gene ID:
84335

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50656-1 targets PRAS40 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRAS40 Monoclonal antibody, PBS Only (Capture) 60477-1-PBS (1F2A9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: PRAS40 Monoclonal antibody, PBS Only (Detector) 60478-1-PBS (1F11C10). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50656-1, PRAS40 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60477-1-PBS. Detection antibody: 60478-1-PBS. Standard: Ag15355. Range: 0.391-100 ng/mL.