

For Research Use Only

KLK3/PSA Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50661-1

Capture Antibody Information

Catalog Number:
60481-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F7C8

Reactivity:
human

GenBank:
NM_001648.2

Conjugate:
Unconjugated

Full name:
kallikrein-related peptidase 3

Gene ID:
354

Detection Antibody Information

Catalog Number:
60481-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3H4B3

Reactivity:
human

GenBank:
NM_001648.2

Conjugate:
Unconjugated

Full name:
kallikrein-related peptidase 3

Gene ID:
354

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-3.125 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50661-1 targets KLK3/PSA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KLK3/PSA Monoclonal antibody, PBS Only (Capture) 60481-1-PBS (1F7C8). 100 μ g. Concentration 1 mg/ml.

Detection antibody: KLK3/PSA Monoclonal antibody, PBS Only (Detector) 60481-2-PBS (3H4B3). 100 μ g. Concentration 1 mg/ml.

Alternative KLK3/PSA matched antibody pairs: MP00975-1, MP00975-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50661-1, KLK3/PSA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60481-1-PBS. Detection antibody: 60481-2-PBS. Standard:Eg0869. Range: 0.098-3.125 ng/mL