

For Research Use Only

PREB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50675-2

Capture Antibody Information

Catalog Number: 66981-2-PBS	Clone ID: 2G11A10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: prolactin regulatory element binding
Isotype: IgG1	GenBank: BC002765	Gene ID: 10113
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28549	

Detection Antibody Information

Catalog Number: 66981-4-PBS	Clone ID: 4B12A8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: prolactin regulatory element binding
Isotype: IgG1	GenBank: BC002765	Gene ID: 10113
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28549	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50675-2 targets PREB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PREB Monoclonal antibody, PBS Only (Capture) 66981-2-PBS (2G11A10). 100 μ g. Concentration 1 mg/ml.

Detection antibody: PREB Monoclonal antibody, PBS Only (Detector) 66981-4-PBS (4B12A8). 100 μ g. Concentration 1 mg/ml.

Alternative PREB matched antibody pairs: MP50675-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50675-2, PREB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66981-2-PBS. Detection antibody: 66981-4-PBS. Standard: Ag28549. Range: 1.563-100 ng/mL.