

For Research Use Only

BCS1L Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50678-1

Capture Antibody Information

Catalog Number:
60212-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3G6H1
Reactivity:
human, mouse
GenBank:
BC007500
Immunogen Catalog Number:
Ag18280

Conjugate:
Unconjugated
Full name:
BCS1-like (yeast)
Gene ID:
617

Detection Antibody Information

Catalog Number:
60212-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1F5H5
Reactivity:
human
GenBank:
BC007500
Immunogen Catalog Number:
Ag18280

Conjugate:
Unconjugated
Full name:
BCS1-like (yeast)
Gene ID:
617

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50678-1 targets BCS1L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: BCS1L Monoclonal antibody, PBS Only (Capture) 60212-1-PBS (3G6H1). 100 µg. Concentration 1 mg/mL.
Detection antibody: BCS1L Monoclonal antibody, PBS Only (Detector) 60212-2-PBS (1F5H5). 100 µg. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50678-1, BCS1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60212-1-PBS. Detection antibody: 60212-2-PBS. Standard: Ag18280. Range: 1.563-200 ng/mL.