

For Research Use Only

FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50689-3

Capture Antibody Information

Catalog Number:
60493-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3G10B1

Reactivity:
human

GenBank:
NM_005860.3

Conjugate:
Unconjugated

Full name:
follistatin-like 3 (secreted glycoprotein)

Gene ID:
10272

Detection Antibody Information

Catalog Number:
60493-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3H2F11

Reactivity:
human

GenBank:
NM_005860.3

Conjugate:
Unconjugated

Full name:
follistatin-like 3 (secreted glycoprotein)

Gene ID:
10272

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50689-3 targets FSTL3/FLRG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FSTL3/FLRG Monoclonal antibody, PBS Only (Capture) 60493-3-PBS (3G10B1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: FSTL3/FLRG Monoclonal antibody, PBS Only (Detector) 60493-5-PBS (3H2F11). 100 μ g. Concentration 1 mg/ml.

Alternative FSTL3/FLRG matched antibody pairs: MP00667-1, MP50689-1, MP50689-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50689-3, FSTL3/FLRG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60493-3-PBS. Detection antibody: 60493-5-PBS. Standard:Eg0949. Range: 0.781-100 ng/mL