

# NINJ2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50693-2

## Capture Antibody Information

**Catalog Number:**  
66202-4-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
5E8C2

**Reactivity:**  
human

**GenBank:**  
BC057766

**Immunogen Catalog Number:**  
Ag5343

**Conjugate:**  
Unconjugated

**Full name:**  
ninjurin 2

**Gene ID:**  
4815

## Detection Antibody Information

**Catalog Number:**  
66202-3-PBS

**Host:**  
Mouse

**Isotype:**  
IgG2b

**Purification Method:**  
Protein A purification

**Clone ID:**  
1D1B10

**Reactivity:**  
human

**GenBank:**  
BC057766

**Immunogen Catalog Number:**  
Ag5343

**Conjugate:**  
Unconjugated

**Full name:**  
ninjurin 2

**Gene ID:**  
4815

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.391-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50693-2 targets NINJ2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NINJ2 Monoclonal antibody, PBS Only (Capture) 66202-4-PBS (5E8C2). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: NINJ2 Monoclonal antibody, PBS Only (Detector) 66202-3-PBS (1D1B10). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative NINJ2 matched antibody pairs: MP50693-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50693-2, NINJ2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66202-4-PBS. Detection antibody: 66202-3-PBS. Standard: Ag5343. Range: 0.391-100 ng/mL.