

# SH3BP5 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50694-3

## Capture Antibody Information

**Catalog Number:**  
67841-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1E6E5  
**Reactivity:**  
human  
**GenBank:**  
BC010123  
**Immunogen Catalog Number:**  
Ag30252

**Conjugate:**  
Unconjugated  
**Full name:**  
SH3-domain binding protein 5 (BTK-associated)  
**Gene ID:**  
9467

## Detection Antibody Information

**Catalog Number:**  
67841-5-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
3D8G8  
**Reactivity:**  
human  
**GenBank:**  
BC010123  
**Immunogen Catalog Number:**  
Ag30252

**Conjugate:**  
Unconjugated  
**Full name:**  
SH3-domain binding protein 5 (BTK-associated)  
**Gene ID:**  
9467

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.391-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.


## Product Information

MP50694-3 targets SH3BP5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: SH3BP5 Monoclonal antibody, PBS Only (Capture) 67841-2-PBS (1E6E5). 100 µg. Concentration 1 mg/mL.  
Detection antibody: SH3BP5 Monoclonal antibody, PBS Only (Detector) 67841-5-PBS (3D8G8). 100 µg. Concentration 1 mg/mL.  
Alternative SH3BP5 matched antibody pairs: MP50694-1, MP50694-2  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50694-3, SH3BP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67841-2-PBS. Detection antibody: 67841-5-PBS. Standard: Ag30252. Range: 0.391-50 ng/mL