

For Research Use Only

# NME3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50704-2

## Capture Antibody Information

<b>Catalog Number:</b> 60505-3-PBS	<b>Clone ID:</b> 1B11E9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> non-metastatic cells 3, protein expressed in
<b>Isotype:</b> IgG2a	<b>GenBank:</b> BC000250	<b>Gene ID:</b> 4832
<b>Purification Method:</b> Protein A Magarose purification	<b>Immunogen Catalog Number:</b> Ag7298	

## Detection Antibody Information

<b>Catalog Number:</b> 60505-4-PBS	<b>Clone ID:</b> 1A2B7	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> non-metastatic cells 3, protein expressed in
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC000250	<b>Gene ID:</b> 4832
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag7298	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.391-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50704-2 targets NME3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NME3 Monoclonal antibody, PBS Only (Capture) 60505-3-PBS (1B11E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: NME3 Monoclonal antibody, PBS Only (Detector) 60505-4-PBS (1A2B7). 100 µg. Concentration 1 mg/ml.

Alternative NME3 matched antibody pairs: MP01200-1, MP01200-2, MP01200-3, MP50704-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50704-2, NME3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60505-3-PBS. Detection antibody: 60505-4-PBS. Standard: Ag7298. Range: 0.391-100 ng/mL.