

For Research Use Only

# FLAD1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50710-1

## Capture Antibody Information

Catalog Number:  
68491-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1C11A10

Reactivity:  
human

GenBank:  
BC011378

Immunogen Catalog Number:  
Ag6845

Conjugate:  
Unconjugated

Full name:  
FAD1 flavin adenine dinucleotide synthetase homolog (*S. cerevisiae*)

Gene ID:  
80308

## Detection Antibody Information

Catalog Number:  
68491-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1D6C2

Reactivity:  
human

GenBank:  
BC011378

Immunogen Catalog Number:  
Ag6845

Conjugate:  
Unconjugated

Full name:  
FAD1 flavin adenine dinucleotide synthetase homolog (*S. cerevisiae*)

Gene ID:  
80308

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50710-1 targets FLAD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FLAD1 Monoclonal antibody, PBS Only (Capture) 68491-2-PBS (1C11A10). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: FLAD1 Monoclonal antibody, PBS Only (Detector) 68491-3-PBS (1D6C2). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50710-1, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-2-PBS. Detection antibody: 68491-3-PBS. Standard: Ag6845. Range: 0.098-100 ng/mL.