

For Research Use Only

FLAD1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50710-3

Capture Antibody Information

Catalog Number:
68491-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F3C6

Reactivity:
human

GenBank:
BC011378

Immunogen Catalog Number:
Ag6845

Conjugate:
Unconjugated

Full name:
FAD1 flavin adenine dinucleotide synthetase homolog (*S. cerevisiae*)

Gene ID:
80308

Detection Antibody Information

Catalog Number:
68491-6-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F10D10

Reactivity:
human

GenBank:
BC011378

Immunogen Catalog Number:
Ag6845

Conjugate:
Unconjugated

Full name:
FAD1 flavin adenine dinucleotide synthetase homolog (*S. cerevisiae*)

Gene ID:
80308

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50710-3 targets FLAD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FLAD1 Monoclonal antibody, PBS Only (Capture) 68491-5-PBS (1F3C6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: FLAD1 Monoclonal antibody, PBS Only (Detector) 68491-6-PBS (1F10D10). 100 μ g. Concentration 1 mg/mL.

Alternative FLAD1 matched antibody pairs: MP50710-1, MP50710-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50710-3, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-5-PBS. Detection antibody: 68491-6-PBS. Standard: Ag6845. Range: 0.781-100 ng/mL.