

For Research Use Only

C1QBP Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50719-1

Capture Antibody Information

Catalog Number:
68084-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
1A8B2
Reactivity:
human
GenBank:
BC013731
Immunogen Catalog Number:
Ag19773

Conjugate:
Unconjugated
Full name:
complement component 1, q subcomponent binding protein
Gene ID:
708

Detection Antibody Information

Catalog Number:
68084-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3A9B2
Reactivity:
human
GenBank:
BC013731
Immunogen Catalog Number:
Ag19773

Conjugate:
Unconjugated
Full name:
complement component 1, q subcomponent binding protein
Gene ID:
708

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50719-1 targets C1QBP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C1QBP Monoclonal antibody, PBS Only (Capture) 68084-2-PBS (1A8B2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: C1QBP Monoclonal antibody, PBS Only (Detector) 68084-3-PBS (3A9B2). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50719-1, C1QBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68084-2-PBS. Detection antibody: 68084-3-PBS. Standard:Ag19773. Range: 0.391-100 ng/mL.